

Is Your Vaccine Management Adequate?



24/7 Temperature Monitoring Protection

Completes Vaccine Temperature Audit Trail

Immediate Problem Alert

Mobile SMS Text and E-mail Alarms

Long Term Temperature Data Traceability

Empowers Existing Fridges

Fridge Door Access Alert and Logs

Standard Spreadsheet Support For Graphs

Up To 16 Fridges Monitored

Trusted Vaccine Storage

Keeping your valuable biologicals safe...

Thanks to NetCommander's 24/7
Fridge Monitoring and Logging Solution
You can be sure your vaccine storage is trustworthy

The NetCommander Vaccine Monitoring Solution

NetCommander can be added to any dispensary fridge to upgrade your vaccine storage and monitoring capabilities. NetCommander is a cost effective web accessible device for monitoring and logging the vaccine storage conditions within dispensary fridges, such as temperature monitoring and ensuring door closure. NetCommander can empower any Bloodbank, Vaccine storage facility, pharmacy fridge, surgery fridge, pathology lab, specimen storage or dispensary fridge.



Primarily aimed at GP surgeries and dispensaries, NetCommander detects changes in storage conditions before vaccines are damaged. The temperature of up to 16 dispensary fridges, or 8 fridges if door monitors are also included, can be monitored and logged simultaneously.

Autonomous 24hr monitoring ensures that your biologicals are maintained at the correct storage temperature even during the night and weekends. User defined logging intervals and 24hr logging provides a comprehensive temperature and status log for all fridges under NetCommander supervision. Periodic log downloads allow easy validation and auditing of storage facilities, allowing you to complete the audit chain from

manufacturer to patient. All logs generated by NetCommander can be imported in to standard spreadsheet applications for graphing and archiving.

Users can access real time temperatures and status logs via their browser, or receive defined alert messages via Windows popup's, e-mail or mobile SMS text messages. For example, if any one of the dispensary fridge doors are left open for more than two minutes an alert can be sent to a designated person. Cameras can be added to the NetCommander solution to add security to restricted access fridges. NetCommander helps you deliver the quality and trust your patients expect.

Technical Specification

Features

- 24/7 monitoring for up to 16 dispensary fridges
- Log temperature every hour
- Ethernet attached standalone monitor
- Wi-Fi option
- No client software required
- Powerful web based interface
- Real time information and status
- Windows, email and SMS text alerts
- On board data storage
- Graphing
- Camera option

Sensors

Temperature	-10-50°C \pm 0.7°C
Humidity	10-90% \pm 6%
Door	

NetCommander General

Operating Temp	0 to 40 °C
Humidity	20-80% non-condensing
Power Supply	240v ac
Dimensions	164mm x 104mm x 55mm
Weight	500g approx

CE approved



CubeRoot Ltd
18 Melton Street
Leicester
LE1 3NB UK

Tel: +44 (0)116 2627370
Fax: +44 (0)116 2628090

www.netcommander.biz

All trademarks acknowledged. Specification subject to change without notice.